




Page
552-564

▪ Suction Instruments



XY-001 - 006

BIEMER

XY-001 4 mm

XY-002 6 mm

XY-003 8 mm

XY-004 10 mm

XY-005 12 mm

XY-006 14 mm

Curettes for aspiration lipectomy, graduated, slightly curved

XY-007 - 008

BIEMER

XY-007 4 mm

XY-008 6 mm

Curettes for aspiration lipectomy, with three apertures, slightly curved malleable

XY-009 - 011

BIEMER

XY-009 8 mm

XY-010 10 mm

XY-011 12 mm

Curettes for aspiration lipectomy, with three apertures, slightly curved malleable

XY-012

For aspirating fluid from the excretory ducts of large glands and for keeping the operating area clear 170 mm

XY-013

Trocar For puncturing 150 mm

XY-014

For clearing large quantities of fluid 190 mm




XY-015
Trocar
For puncturing
170 mm

XY-016
Hollow needle,
sharp
150 mm

XY-017
double curved
For general
surgery
230 mm

XY-018
For tonsillectomy,
oral and laryngeal
operations
190 mm

XY-019
For tonsillectomy,
oral and laryngeal
operations
220 mm





XY-026 8 CH
 XY-027 10 CH
 XY-028 12 CH
 XY-029 14 CH
 XY-030 16 CH
 XY-031 18 CH

REDON
 XY-026 - 031
 threaded end
 19.5 cm



XY-032 8 CH
 XY-033 10 CH
 XY-034 12 CH
 XY-035 14 CH
 XY-036 16 CH
 XY-037 18 CH

REDON
 XY-032 - 037
 threaded end
 14 cm



XY-038 8 CH
 XY-039 10 CH
 XY-040 12 CH
 XY-041 14 CH
 XY-042 16 CH
 XY-043 18 CH

REDON
 XY-038 - 043
 threaded end
 19.5 cm



XY-044 8 CH
 XY-045 10 CH
 XY-046 12 CH
 XY-047 14 CH
 XY-048 16 CH
 XY-049 18 CH

REDON
 XY-044 - 049
 threaded end
 19.5 cm



POOLE
XY-050
24 cm
Ø10 mm



POOLE
XY-051
25 cm
Ø8 mm




COLLEY
XY-052
with perforated screw tip
31 cm
Ø10 mm




COLLEY
XY-053
with perforated screw tip
31 cm
Ø7 mm



COLLEY
XY-054
36 cm
Ø8 mm



COLLEY
XY-055
with perforated screw tip
36 cm
Ø8 mm





1/2

**YANKAUER
XY-069**
complete, with suction
tip and connector
31 cm



1/2

**YANKAUER
XY-070**
complete, with handle, 4 tube,
suction tip and tube connector
31 cm



1/2

**Fig. 0
XY-073**
10 cm



**Fig. 1
XY-074**
10 cm



**Fig. 2
XY-075**
10 cm



**Fig. 3
XY-076**
complete, with
suction tip
13.5 cm



1/2

XY-071
Handle
15 cm



1/2

XY-072
Olive
2.5 cm

WULLSTEIN

XY-077
 Suction tube, convical,silver,oliver
 with LUER hub
 16 cm

Ø1,75 mm

Ø3 mm

 1/2
 Ø8 mm

WULLSTEIN

XY-078 1.5 mm**XY-079** 2 mm
 Suction tube, cylindrical,silver,with
 stylet with LUER hub
 16 cm

 1/2
 Ø8 mm

Usable length

YASARGIL usable length diam.

XY-080	13 cm	1.5 mm
XY-081	15 cm	1.5 mm
XY-082	18 cm	1.5 mm
XY-083	22 cm	1.5 mm
XY-084	13 cm	2 mm
XY-085	15 cm	2 mm
XY-086	18 cm	2 mm
XY-087	22 cm	2 mm

YASARGIL usable length diam.

XY-088	13 cm	2.5 mm
XY-089	15 cm	2.5 mm
XY-090	18 cm	2.5 mm
XY-091	22 cm	2.5 mm
XY-092	13 cm	3.5 mm
XY-093	15 cm	3.5 mm
XY-094	18 cm	3.5 mm
XY-095	22 cm	3.5 mm

 Ø6.5-9 mm
 with LUER hub

LANDOLT

XY-096 - 098

15 cm

diam Ø mm

XY-096 8 mm / curved**XY-097** 11mm / curved**XY-098** 15mm / curved

Ø3.5 mm





Ø2.0 mm

XY-099 5 cm**XY-100** 35 cm
 with rubber tip,
 with LUER hub

Ø9 mm

XY-101 2 mm diam**XY-102** 2.5 mm diam**XY-103** 3 mm diam
 with finger cut-gff,
 thumb piece 90° to
 curved, with LUER hub
 25 cm

 Ø6.9 mm
 1/2





1/2

**WULLSTEIN
XY-112 - 114**
14 cm
handle with tube, only



**WULLSTEIN
XY-112**
4.5 cm
Mirco Suction Cannula
Complete, with handle and 3 tips
Ø 0.4 mm



**WULLSTEIN
XY-113**
4.5 cm
Mirco Suction Cannula
Complete, with handle and 3 tips
Ø 0.6 mm



**WULLSTEIN
XY-114**
4.5 cm
Mirco Suction Cannula
Complete, with handle and 3 tips
Ø 0.8 mm



XY-115
working length 75 mm
Ø 1 mm



XY-116
working length 75 mm
Ø 1.5 mm



XY-117
working length 75 mm
Ø 2 mm



1/2

**BARON-SCHUKNECHT
XY-115 - 117**
with finer cut-off



1/2

**ZOELLNER
XY-118 - 121**
Suction Cannula


Suction Tip

XY-118 Ø 0.5 mm
XY-119 Ø 0.7 mm
XY-120 Ø 0.9 mm
XY-121 Ø 1.2 mm



LUER
XY-122 - 124
Conical

1/2



XY-122
16.5 cm
Ø 3 mm
working length
110 mm



XY-123
16.5 cm
Ø 3.5 mm
working length
110 mm



XY-124
18 cm
Ø 4.5 mm
working length
130 mm



YASARGIL
XY-125 - 128
working length
15 cm

1/2




XY-125
17 cm
Ø 1.5 mm




XY-126
17 cm
Ø 2 mm



XY-127
17 cm
Ø 2.5 mm



XY-128
17 cm
Ø 3.5 mm





Ø 6-9 mm
XY-135 Ø 1,5 mm
XY-136 Ø 2,0 mm
XY-137 Ø 2,5 mm
XY-138 Ø 3,0 mm

1/2
XY-135 - 138
PLESTER
with finger cut-off and
stylet, thumb piece
parallel to curved, with
LUER hub
195 mm



Ø 6-9 mm
XY-139 2,0 mm
XY-140 2,5 mm
XY-141 3,0 mm
XY-142 4,0 mm
XY-143 5,0 mm

1/2
XY-139 - 143
FERGUSON
with finger cut-off and
stylet, thumb piece
90° to curved, with
LUER hub
180 mm



Ø 6-9 mm
XY-144 Ø 4,0 mm
FERGUSON
insulated, with finger cut-off
and stylet, thumb piece 90°
to curved, with LUER hub
180 mm



Ø 6-9 mm
XY-145 Ø 3,0 mm
XY-146 Ø 5,0 mm
FRAZIER
insulated, with finger
cut-off and stylet, thumb
piece 90° to curved,
with LUER hub
180 mm



1/2
XY-147
5,0 mm diam., Tip
6,5 mm wide, Nasal
suction tube, with
LUER hub
180 mm



1/2
XY-148 Ø 4,5 mm
XY-149 Ø 5,5 mm
Nasal suction
tubes
175 mm



1/2
XY-150
GORNEY
Suction dissector
180 mm



Ø 6-10 mm
1/2
XY-151 4 mm wide
XY-152 5 mm wide
Suction raspatories,
sharp, bayonet-shaped,
with stylet
200 mm



Ø 9 mm
1/2
XY-153
LOPEZ-REINKE
Tonsil dissector with
suction device
235 mm



XY-154
SPACKMAN CANNULA
Surgical and Latexfree
Cone the depth of the uterus
and to allow
for gentle dilation of the cervix
Length: 360 mm