


Page 125-140




▪ Orthodontic Pliers



INSERTED PLIERS



NON-INSERTED PLIERS




23-1001-000
 Young Plier

23-1002-000
 Laboratory Plier
 Smooth

23-1003-000
 Robinson Plier

23-1004-000
 Reynolds Plier

23-1005-000
 How Plier

23-1006-000
 Goslee Plier

23-1007-000
 Laboratory Plier
 Serrated

23-1008-000
 Pesso Plier

23-1009-000
 Johnson Plier

23-1010-000
 Three Jaw Plier-Large

23-1011-000
 Gordon Plier

23-1012-000
 Three Jaw Plier-Large
 Curved


Curved


23-1013-000
 Concaved Clasp
 Bending Pliers, Straight


Straight

BOX SELECTION



Box Lock=B/L




Screw Lock Box=S



Allen Wrench Box=W

Clasp Foring Pliers




22-1014-000
Type-I



22-1015-000
Type-II



22-1016-002
Type # 2



22-1017-000
Type-III



22-1018-000
Type-IV



22-1019-000
BAR BENDERS
Type # 1



22-1020-000
BAR BENDERS

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/Allen Wrench Box=W/L






23-1021-TC
Distal End Safety Cutter
-Slim (Long Handle)
Maximum Cutting Capacity
■ Rectangular Wires up to .022 x .028"
● Round Wires up to .020"
♥ Brided Wires up to .016"




23-1022-TC
Distal End Safety Cutter
-Slim
Maximum Cutting Capacity
■ Rectangular Wires up to .022 x .028"
● Round Wires up to .020"
♥ Brided Wires up to .016"



By regular distal end
safety cutter

By distal end safety
cutter flush cut type



23-1024-TC
Distal End Safety Cutter
-Slim (Flush Cut Type)
Maximum Cutting Capacity
■ Rectangular Wires up to .022 x .028"
● Round Wires up to .020"
♥ Brided Wires up to .021 x .025



23-1023-TC
Distal End Safety Cutter
Maximum Cutting Capacity
■ Rectangular Wires up to .022 x .028"
● Round Wires up to .020"
♥ Brided Wires up to .016"



23-1025-TC
Pin & Ligature Cutter
Maximum Cutting Capacity
Ligature Wires up to .016"
Elastic Ligatures



23-1026-TC
Distal End Cutter
-Non Hold
Maximum Cutting Capacity
■ Rectangular Wires up to
.022 x .028"



23-1027-TC
Pin & Ligature Cutter-Miniature
Maximum Cutting Capacity
Ligature Wires up to .016"
Elastic Ligatures



23-1028-TC
Pin & Ligature Cutter
-Slim (Long Handle)
Maximum Cutting Capacity
Ligature Wires up to .015"
Elastic Ligatures

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/Allen Wrench Box=W/L





23-1029-TC
Pin & Ligature Cutter
Maximum Cutting Capacity
Soft Wires up to .015"



23-1030-TC
Pin & Ligature Cutter-Miniature
Maximum Cutting Capacity
Soft Wires up to .015"



23-1031-TC
Pin & Ligature Cutter-Slim
Maximum Cutting Capacity
Soft Wires up to .015"




23-1032-TC
Pin & Ligature Cutter-45° Angel
Maximum Cutting Capacity
Soft Wires up to .015"




23-1033-TC
Hard Wire Cutters
Maximum Cutting Capacity
Hard Wires up to .020" Round



23-1034-TC
Pin & Ligature Cutters
Maximum Cutting Capacity
Ligature Wires up to .016"
Elastic Ligature



23-1035-000
Power Cutters
Maximum Cutting Capacity
Hard Wires & Soft Wires
up to 1.5mm



23-1036-TC
Wire Hold Cutters
**Maximum Cutting Capacity
& Holding Capacity**
Hard Wires of 0.7mm to 1.2mm



23-1037-000
Flat Cutters
Maximum Cutting Capacity
Cut at Top ~ 1.00mm
Cut at Side ~ 1.2mm

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/Allen



Wrench Box=W/L



23-1038-000
How Plier-Titanium Alloy
With Special coated beaks



23-1039-000
How Plier
Tips are 1/8" (3.2mm)
diameter



23-1040-000
How Plier
Tips are 1/8" (3.2mm)
diameter



23-1041-000
How Plier-Small
Tips are 3/32" (2.4mm)
diameter



23-1042-000
How Plier-Offset
Tips are 1/8" (3.2mm)
diameter



23-1043-000
Weingart Utility Plier
-Titanium Alloy
with special coated beaks



23-1044-000
Weingart Utility Plier
New Type of the Weingart Utility Plier
designed for a more precise operation.



23-1045-TC
Weingart Utility Plier



23-1046-000
Weingart Utility Plier



23-1047-TC
Weingart Utility Plier



23-1048-TC
Bird Beak Plier with Cutter
**Maximum Bending &
Cutting Capacity**
Wires upto .032"



23-1049-TC
Bird Beak Plier
**Maximum Bending
Capacity**
Wires upto .032"



23-1050-000
Bird Beak Plier
**Maximum Bending
Capacity**
Wires upto .030"

BOX SELECTION




Box Lock=B/L



Screw Lock Box=S/L Allen Wrench Box=W/L





23-1051-TC
Light Wire Plier with Cutter
-Grooved Square Tip
**Maximum Bending &
Cutting Capacity**
Wires up to .020"



23-1052-TC
Light Wire Plier with Cutter
-Grooved Square Tip
Maximum Bending Capacity
Wires up to .020"



23-1053-000
Light Wire Plier with Cutter
-Grooved Square Tip
Maximum Bending Capacity
Wires up to .018"



23-1054-TC
Light Wire Plier with Cutter
**Maximum Bending &
Cutting Capacity**
Wires up to .020"



23-1055-TC
Light Wire Plier
**Maximum Bending
Capacity**
Wires up to .020"



23-1056-000
Light Wire Plier
**Maximum Bending
Capacity**
Wires up to .018"



23-1057-000
Begg Light Wire Plier
**Maximum Bending
Capacity**
Wires up to .016"



23-1058-000
Young Plier
**Maximum Bending
Capacity**
Wires up to .040"



23-1059-000
Jarabak Plier
**Maximum Bending
Capacity**
Wires up to .020"



23-1060-000
Tweed Loop Forming Plier
**Maximum Bending
Capacity**
Wires up to .022"



23-1061-000
Omega Loop Forming Plier
**Maximum Bending
Capacity**
Wires up to .022"



23-1062-000
Loop Forming Plier
**Maximum Bending
Capacity**
Wires up to .022"

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S



Allen Wrench Box=W/L



23-1063-TC
 Quad Helix Bending Plier
Maximum Bending Capacity
 Wires up to .036"



23-1064-TC
 Tweed Loop Forming Plier
Maximum Cutting Capacity
 Rectangular Wires up to
 .022 x .028"



23-1065-TC
 Kim Plier with Cutter
Maximum Cutting Capacity
 Rectangular Wires up to
 .016 x .022"



23-1066-000
 Nancy Closing Loop Plier
Maximum Bending Capacity
 Wires up to .022"



23-1067-000
 Tweed Arch Bending Plier
Maximum Cutting Capacity
 Rectangular Wires up to
 .022 x .028"



23-1068-TC
 Tweed Arch Bending Plier
Maximum Cutting Capacity
 Rectangular Wires up to
 .022 x .028"



23-1069-TC
 Three Jaw Plier
Maximum Bending Capacity
 Wires up to .028"



23-1070-000
 Three Jaw Plier
Maximum Bending Capacity
 Wires up to .028"



23-1071-000
 Heavy Three Jaw Plier
Maximum Bending Capacity
 Wires up to .036"



23-1072-000
 Adams Plier
Maximum Bending Capacity
 Wires up to .030"



23-1073-000
 Adams Plier
Maximum Bending Capacity
 Wires up to .040"



23-1074-000
 Stop Plier
Maximum Bending Capacity
 Wires up to .025"



23-1075-TC
 Heavy Wire Bending Plier
Maximum Bending Capacity
 Wires up to .047" (1.2mm)



23-1076-000
 Adams Plier
Maximum Bending Capacity
 Wires up to .028" (0.7mm)

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S



Allen Wrench Box=W/L



23-1077-000
Head-Gear Plier
Maximum Bending Capacity
Wires up to .051"(1.3mm)



23-1078-000
Head-Gear Plier
Maximum Bending Capacity
Wires up to .051"(1.3mm)



23-1079-TC
Angel Three Jaw-Left




23-1080-TC
Angel Three Jaw-Right




23-1081-TC
Arch Forming Plier
-Grooved



23-1082-TC
Arch Forming Plier
-Non Grooved



23-1083-TC
Ligature Forming Plier



23-1087-000
Bracket Placing
Tweezers



23-1084-000
Hollow Chop Plier
Maximum Bending Capacity
Wires up to .022"



23-1085-000
Tube Crimping Plier



23-1086-000
Hook Crimping Plier
For use in quick and precise
placement of the hook on the archwire.



23-1088-000
Boon Bracket Positioning
Gauge

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/Allen Wrench Box W/L





23-1089-001
Bracket Holding
Tweezers,
Scraper End



23-1090-002
Bracket Holding
Tweezers,
Round



23-1091-003
Bracket Holding
Tweezers,
Round End



23-1092-004
Bracket Holding
Tweezers,
with Bracket Aligner



23-1093-005
Bracket Holding
Tweezers,
Curved




23-1094-006
Bracket Holding
Tweezers,
45 degree




23-1095-007
Bracket Holding
Tweezers,
90 degree



23-1096-008
Buccal Tube
Tweezers




23-1097-TC
Direct Bond Attachment
Removing Plier-Angled



23-1098-TC
Direct Bond Bracket Attachment
Removing Plier-Straight



23-1099-TC
Direct Bond Bracket Attachment
Removing Plier-Straight



23-1100-TC
Direct Bond Attachment
Removing Plier-Angled



23-1101-TC
Bond Removing Plier



23-1102-000
DBS Bracket Removing Plier



23-1103-000
Ligature Typing Plier
-Coon Type(Safety Type)



23-1104-000
Ligature Typing Plier



23-1105-000
Mathieu Type

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/L



Allen Wrench Box=W/L




23-1106-000
Elastic Placing Plier



23-1107-000
Ligating Plier



23-1108-000
Ligature Utility Plier



23-1109-000
Elastic Placing Plier
Hook Tip



23-1111-000
Mosquito
Forceps
Curved

23-1110-000
Mosquito Forceps-Straight



23-1112-000
Ligature Director



23-1113-000
Ligature & Director



23-1114-000
Twister



23-1115-000
Ligature Director-Single End



23-1116-001
Distal End Bender/Twister, .22" (.56mm)

23-1117-002
.30" (.76mm)



23-1118-000
Band Contouring Plier



23-1119-000
DE LA Rosa Plier



23-1120-000
Crown Contouring Plier



23-1121-000
Crown Crimping Plier



23-1122-000
Crown Crimping
Plier-Slim

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/L



Allen Wrench Box=W/L



23-1123-TC
Posterior Band Removing
Plier



23-1124-TC
Posterior Band Removing
Plier



23-1125-TC
Posterior Band Removing
Plier-Long



23-1126-000
Posterior Band Removing
Plier



23-1127-TC
Posterior Band Removing
Plier-Long



23-1128-TC
Posterior Band Removing
Plier



23-1129-000
Anterior Band Removing
Plier



23-1130-000
Ki-Type Bar
Bending Pliers
Type # 2



23-1131-000
Big Young Plier
Bending Plier



23-1132-000
Ki-Type Bar Bending Pliers
Type # 1

BOX SELECTION



Box Lock=B/L



Screw Lock Box=S/L



Allen Wrench Box=W/L



23-1133-000
Direct Bond Bracket



23-1134-000
Wire Bending Plier



23-1135-TC
Wire Twister 6"
with Carb-N-Sert Jaws



23-1136-TC
Wire Twister 8"
with Carb-N-Sert Jaws